

Primary cell line antibiotic sensitivities ($\mu\text{g/ml}$)

Primary GBM cell line	Puromycin	Hygromycin B	G418	Blasticidin S
BAH1	0.5 (5)	100	200	1
FPW1	0.5	200 (9)	300 (9)	4
HW1	0.4	200	300	2
JK2	0.4	100	300	2
MMK1	0.5	100	200	2
MN1	0.4 (5)	200	200	2
PB1	0.2 (5)	200	200	1
RK11	0.5	200 (10)	200 (10)	4
RN1	0.5	100	200	1
SB2b	0.5 (5)	200	400	1
SJH1	0.5	200	300	1
WK1	0.5	200 (9)	400	4

Notes: 7 days unless specified ().